



Genlantis, a division of Gene Therapy Systems, Inc.
10190 Telesis Court, San Diego, California 92121
Toll Free: 888-428-0558. **Telephone:** 858-457-1919
Fax: (858) 642-0937
Email: orders@genlantis.com
URL: www.genlantis.com

Certificate of Analysis

PrimaPure™ Bovine Smooth Muscle Cell Growth Medium was designed to promote attachment, spreading and proliferation of bovine smooth muscle cells in culture. This medium is fully supplemented with fetal bovine serum, growth factors and antibiotics. Bovine Smooth Muscle Cell Growth Medium is sterile filtered, bicarbonate buffered, and ready to be used for culturing bovine smooth muscle cells.

Product Name: PrimaPure™ Bovine Smooth Muscle Cell Growth Medium
Quantity: 500 ml
Cat. No.: PMB311500

Quality Control:

Bioassay: Each lot is tested for its ability to promote plating and proliferation of bovine smooth muscle cells in culture.
Sterility Test: Negative for bacteria, yeast and fungi.

Storage and Stability: Avoid exposure to light. Stable for 2 months at 4°C.