



Genlantis, a division of Gene Therapy Systems, Inc.
10190 Telesis Court, San Diego, California 92121
Toll Free: 888-428-0558. **Telephone:** 858-457-1919
Fax: (858) 642-0937
Email: orders@genlantis.com
URL: www.genlantis.com

Certificate of Analysis

This Human Microvascular Endothelial Cell (HMVEC) Growth Medium is a new and improved version of growth medium for promoting attachment, spreading and proliferation of HMVEC on Attachment Factor Solution-coated tissue culture ware. This medium is sterile filtered and fully supplemented with fetal bovine serum, growth factors, trace elements, antibiotics, and is ready for use. HMVEC Growth Medium does not contain PDGF or VEGF.

Product Name: PrimaPure™ HMVEC Growth Medium
Quantity: 500 ml
Cat. No.: PM112500

Quality Control:

Bioassay: Each lot is tested for its ability to promote the attachment, spreading and proliferation of HMVEC on Attachment Factor Solution coated tissue culture ware.
Sterility Test: Negative for bacteria, yeast and fungi.

Storage and Stability: Avoid exposure to light. Stable for 2 months at 4°C.