



Genlantis, a division of Gene Therapy Systems, Inc.
10190 Telesis Court, San Diego, California 92121
Toll Free: 888-428-0558. **Telephone:** 858-457-1919
Fax: (858) 642-0937
Email: orders@genlantis.com
URL: www.genlantis.com

Certificate of Analysis

PrimaPure™ Smooth Muscle Cell Growth Medium was designed to promote attachment, spreading and proliferation of smooth muscle cells in culture. This medium is fully supplemented with amino acids, vitamins, minerals, organic compounds, inorganic salts, trace elements, fetal bovine serum, growth factors, and antibiotics. PrimaPure™ Smooth Muscle Cell Growth Medium is sterile filtered, bicarbonate buffered, and is ready for use.

Product Name: PrimaPure™ Smooth Muscle Cell Growth Medium
Quantity: 500 ml
Cat. No.: PM311500

Quality Control:

Bioassay: Each lot is tested for its ability to promote plating, spreading and proliferation of smooth muscle cells in culture.
Sterility Test: Negative for bacteria, yeast and fungi.

Storage and Stability: Avoid exposure to light. Stable for 2 months at 4°C.